

testo DiSCmini data conversion tool

This software tool loads and analyzes raw DiSCmini data files which are saved on the SD card.

It allows to view the data and calculate particle number concentration, mean diameter and LDSA. You may also export calculated data for using it in other applications like Microsoft Excel.

Make sure a Java Runtime Environment or Java Development Kit is installed on your PC and just double-click the file testo_DiSCmini_dct.bat

For further information:

- see testo DiSCmini data conversion tool -> menu "Help"
- or contact our technical support: support-nanoparticle@testo.de
(please include screenshots of potential error messages)